



Product description

Order number	87.7708.800
Product description	High-performance RF sealing device, with digital display, for blood bag tubes up to 6.4 mm Ø, with optical tube detection, with touchscreen operation, incl. transport case

Product characteristics

Application	Obtaining blood and blood products
Supply voltage	230 V
Power consumption during operation	200 W
Frequency	50 Hz - 60 Hz

Size

Width of product	7,2 cm
Height of product	24 cm
Length of product	34 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Colour of product grey

Purity & certification

Product category no medical device | IVD
CE certified CE
Batched no

Packaging

Minimum order qty. 1
Piece(s) / inner box 1
Piece(s) / case 1
Piece(s) / pallet 1
Case volume 0,028 cbm
Weight of product 7,147 kg
Weight of case 7,76 kg
EAN of inner box 4038917341608
EAN case 4038917341592

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

