

TS3712-Station

Automatic sealer



Joint R&D Sino-German Technology team

8-year Joint R&D Sino-German team ensure fast and reliable product quality

Modern design

Ergonomic design offers rounded edges for comfortable handling

Light and easy to carry

Portable device for rapid thermal closing of medical PVC conduits and powered by batteries. Sufficient batteries for more than 1500 hot closures of PVC conduits up to 5 mm in diameter

TS3712-Station Automatic Sealing

Best choice of Apheresis center

- High frequency heating method
- Automatic sealing
- Sealing time: 0.6-3sec
- LED Indicator
- Sealing tubes with up to 5mm outer diameter



Technical Information

Device Name	TS3712-Station
Heat Sealing Time	1-3seconds
Dimension	L*W*H:295mm*180mm*130mm
Weight	3.0kg
Voltage	110V-240V
Frequency	50/60Hz